




HOW TO OPERATE



- To switch on your hoist press the grey handset button
- The handset has 2 functions - lift & lower
- Minimal effort from your carer is required to gently guide you along the track



HOW TO CHARGE



- Your hoist charges up when the handset is replaced back into the docking station.
- When placed into the docking station you should hear 3 beeps to confirm the hoist is charging.
- Always keep your handset in the docking station when your hoist is not in use.
- The hoist will automatically switch off when not in use

SAFETY FEATURES



- **EMERGENCY STOP** - in the event of an emergency immediately stop the hoist by pulling the red chord.
- Once the red chord has been pulled, emergency lowering can be activated by gently pulling down on the red chord. The hoist will lower down very slowly.
- If the red chord has been pulled in an emergency or **BY ACCIDENT** the hoist will stop working.
- To re-start the hoist push the yellow button in the hoist pod back in and press any button on the handset. The hoist light should return to green & hoist will start working again.




SAFETY FEATURES

- LOWERING SAFETY FEATURE - if an object is obstructing the carry bar an anti-slack safety device prevents the hoist from lowering. Remove the obstruction & the hoist will work again.



MAINTENANCE

- Wipe clean with a damp cloth & keep rails free from dust. Do not use cleaning solvents.
- SERVICING - your hoist should be serviced every 6 months - look for a white service sticker on the hoist to see when your next service is due.



www.simplymovingandsleeping.co.uk



info@simplymovingandsleeping.co.uk



01422 5812 34 / 01913 941234

TROUBLESHOOTING

GH1 HOIST NOT WORKING

STEP 1

CHECK POWER

Is the lead to the docking station plugged in?

YES

NO - PLUG IN TO RESOLVE

Is there a green light showing on the right side of the docking station?

YES - Is there a light showing on the under side of the hoist?

YES - MOVE TO STEP 2

NO

Press the handset UP or DOWN button - has the light come on now?

YES - MOVE TO STEP 2

NO - MOVE TO STEP 3

STEP 2

CHECK CHARGING

What colour is the charging light on the under side of the hoist ?

GREEN

Replace the handset into the docking station - listen for 3 beeps

MOVE TO STEP 3

AMBER

Hoist needs to charge - leave your handset in the docking station & leave for 20 mins to charge - **ISSUE SHOULD RESOLVE**

****TIP** - always leave your handset plugged into the docking station after using your hoist

STEP 3

CHECK EMERGENCY STOP

On the underside of the hoist can you see the yellow reset button pushed out?

YES - I can see yellow

Push the yellow reset button fully back in

Press the handset UP or DOWN button - **ISSUE SHOULD RESOLVE**

NO

All checks completed & **HOIST STILL NOT WORKING**

CONTACT SIMPLY FOR ADVICE

www.simplymovingandsleeping.co.uk



info@simplymovingandsleeping.co.uk



01422 5812 34 / 01913 941234